

Research Genetic Cancer Centre

RGCC INTERNATIONAL GmbH
HEADQUARTERS
Zug City Centre
2nd Floor
Baarerstrasse 95
Zug 6300
Switzerland / Schweiz / Suisse

Personalised Cancer Testing

- Early detection
- Personalised treatment of cancer
- Chemosensitivity testing

R.G.C.C. International GmbH offers a range of individualised tests from a world-class molecular oncology laboratory in Greece.

These tests can:


- 1) detect early signs of a developing cancer,
- 2) help to monitor existing cancers and
- 3) produce an individual profile of which cancer drugs and which natural substances can be used to achieve the best treatment outcome.

Circulating Tumour Cells

CTCs are cancer cells, which have broken away from the primary tumour and have entered the blood stream where they circulate and have the potential to generate metastatic disease. These cells can be isolated and identified, and there is growing interest in their detection for the following purposes:

- The early detection and diagnosis of new cancers
- Monitoring of existing cancers
- Prognosis providing information about the risk of recurrence of a current or old cancer

R.G.C.C. International GmbH offers general assays for CTCs and for specific cancers, including breast, colorectal, prostate, malignant melanoma and sarcomas.



Chemosensitivity testing

Cancer doctors principally rely on the statistical analysis of large treatment trials, to decide which drugs to use for specific cancers. There is a growing interest in personalized cancer therapy, which involves identifying those treatments, which may work best for an individual's cancer. Chemosensitivity testing is one method of doing this.

Chemosensitivity testing involves testing an individual's cancer cells in the laboratory to see which drugs demonstrate the best response. It therefore provides guidance about which treatments may be best for the individual in clinical practice.

The R.G.C.C. Internation GmbH test is a blood test (or sometimes tissue). Tumor cells are identified and isolated from the sample for the following analysis:

- Viability testing of chemotherapy drugs
- Genetic profiling for guidance about targeted therapies eg monoclonal antibodies
- Viability testing (and identification of mechanisms of action), of natural substances which may be used as part of a complementary treatment strategy.

The results are presented in a written report which your doctor can use to help guide your treatment options and choices.

In addition, R.G.C.C. International GmbH can provide information about how an individual will 'handle' specific chemotherapy agents. Our genetic makeup determines whether we are 'accumulators' or 'rapid metabolizers' of certain drugs. This can play a critical role in determining how effective a specific treatment is likely to be for us, and how significant the side effects will be.

Who are these tests for?

- People who want to actively engage in reducing their risk of developing cancer in the future
- People with an increased risk of cancer e.g. due to family history or lifestyle / environ mental issues who want the opportunity to engage in a screening programme for early detection and diagnosis.
- People with a current diagnosis of cancer who want more information about treatment options for them as an individual – including natural treatments.

Who is R.G.C.C. International GmbH?


R.G.C.C. International GmbH is a world-class laboratory which specialises in medical genetics and in particular cancer genetics, providing a global service. The Director and founder is Dr Ioannis Papasotiriou MD.


Scientific Evidence

- It is well established by literature and validation that the overall response rate from empirical chemotherapy varies from 5 to 7.5% (Royal North Shore Hospital Clin Oncol (R Coll Radiol) 2005 Jun;17(4):294).
- Pharmaceutical industry assesses and validates candidate drugs by utilizing extensive chemosensitivity assays under the term of High Throughput Screening (HTS).
- Cancer is caused by severe damaged genetic material which leads to random genetic instability. So, each malignancy behaves differently to each individual.
- Each person has different genetic fingerprinting from the others. This classifies people to:
 - o Rapid metabolizers (individuals that metabolize fast the drug without receive any benefit from it).
 - o Accumulators (individuals that cannot excrete the byproducts of a drug and that causes severe side effects and toxicities).
 - o Normal metabolizers (these individuals can normally metabolize a drug to its active form and excrete the byproducts normally).

Hence based on the above, each patient and each malignancy has its own identity and behaves individually and differently from person to person.

Therefore personalized treatment becomes essential in order to generate better rates of successful treatments in cancer. (DIAGRAMM: NSCLC of case treated with empirical lines vs personalised options)








Therapeutic Concept

The tumor is consisted from several subgroups of cells (subpopulations) with different features. One of the subset actually drives the progress of the disease, the resistance to therapy and the relapse. This subset is called cancer stem cell like cells or tumor initiating cells. When a patient is treated and the cancer cells may be destroyed, then when the compatible diagnostic does not discover any signal, the cancer cells may consist a population of 10⁹ to 10¹² cells. This limit defines the remission stage of a patient. At that stage only the cancer stem cells may survive and colonize into distant organs and generate metastases in time.

Hence the usage of our test is to detect in which therapeutic approach the cancer cell may respond and also during remission, it is essential to detect, discover and explore the features of the disseminated cancer stem cell like, in order to delay further the risk of relapse and also to generate options to treat even these kinds of cells. The main goal is to discover, analyze and screen the cancer cells in every step of the disease.





Chemosensitivity / Chemoresistance

Oncostat Plus former TUP Plus

Chemosensitivity for Cytotoxic drugs, targeted therapies (MOAB & TKIs) & Natural substances.

This test includes the chemosensitivity / chemoresistance assessment for cytotoxic drugs, monoclonal antibodies, small molecules that inhibit specific targets (TKI etc), the assessment of natural substances and extracts for an anticancer potency.





Samples at: http://www.rgcc-genlab.com/?tests

* CTCs: Circulating Tumor Cells

Oncostat former TUP

Chemosensitivity for Cytotoxic drugs

This test includes the chemosensitivity / chemoresistance assessment for cytotoxic drugs, monoclonal antibodies, small molecules that inhibit specific targets (TKI etc).

Oncostat Extracts

former TUP Extracts

Chemosensitivity only for Natural Substances

This test includes only the assessment of natural substances and plant extracts for potency against cancer patients derived CTCs.

Screening & Follow Up Tools

Oncocount

CTC Number Only

This test has as an outcome only the enumeration of CTCs, it does not include any other information concerning CTCs.

Oncotrail

CTC Number & Immunophenotype for specifc types of malignancies

This test is a tailor made test for specific type of malignancies such as breast cancer (Oncotrail for breast), prostate cancer (Oncotrail for prostate), Sarcoma etc. This test includes only relevant markers for a specific type of malignancy which makes the test a good tool for follow up control.

Oncotrace

CTC Number & Immunophenotype for all types of malignancies

This test has as an outcome the concentration of CTCs and also the immunophenotype control of this kind of cells.

Immunostat

Profile of humoral and cellular immunity and cachexia

This test uses specific cellular markers and cytokine production to detect the type or types of cells that are responsible for the activation or repression of the immune system of a patient.

Two different assays are used for the result export.

ELISA assay is used for the detection of cytokines produced by specific cells in serum sample.

On the other hand flow cytometry is used for the enumeration of specific markers on peripheral blood cells.

Both assays are then compared to a control sample and the results are represented as the following table:

Immunity activation/suppression

Hematopoietic cells	Cell Type		Markers / cytokines		Control	Results	
CD45 positive cells	B cells	Plasma cells	CD19				个
		Naïve B cells	CD20				\downarrow
		Memory cells	CD28				-
	CD8 (Cytotoxic Lymphocytes, CTI		CD62L				
	T cells	CD4 (T helper cells, $T_{\rm H}$)	<u>Th1 (cell</u> mediated <u>immunity)</u>	IFN-γ (interferon gamma)			
			<u>Th2 (humoral</u> <u>immunity)</u>	IL-4 (Interleukin 4)			
				IL-6 (Interleukin 6)			
				IL-2 (Interleukin 2)			
			Th3 (suppressed immunity)	TNF-α (Tumor Necrosis Factor alpha)		
				T regulatory CD2:	;		
				cells (Tregs) CTLA	4		
		Memory cells	CD28				
	Antigen		CD80				
	Presenting Cells (APCs)	Dendritic cells (DCs)	CD86				

This test will be a useful tool for practicians dealing with patients with acute and chronical immune diseases.

Metastat

Markers on CTCs that point out the potent organ for relapse

From literature numerous markers are associated with specific tumor organ metastases. The Metastat test is a collection of specific markers that are detected on specific organs.

Cells of the primary tumor that have the ability to metastasize, will express specific genes and proteins according to the organ that will metastasize, and that is given by the genetic information.

This test uses patient's blood samples that are analyzed with two different methods for the export of accurate results.

Real Time PCR detects the levels of these specific marker genes and flow cytometry the protein levels produced by these exact genes.

Results are compared to a control sample and the results are given with the following format

Metastasis Markers

Metastasis location	<u>Marker</u>	Sample levels	Normal levels	<u>Results</u>
	TGF-β R2			
General	ITGB-4 R			V
General	ITGB-5 R			
	ITGB-6 R			
Pleura	CCR6			
Pieura	Mesothelin			
Skin	CCR7			
Lung	IGF-R2			
Lung	Phospho-ERK1/2			
	Bone Morphogenetic			
Bone	Protein Receptors			
Done	CXCR4 and RANK			
	BST-2			
	CXCR4			
Liver	TRAIL-R2			
Livei	FAS R			
	HGFR			
	Phospho-STAT-3			
Brain	CX3CR1			
	DSC-2			

Metastat is a tool that will help practicians in the prognosis of metastases trends in cancer patients and guide them in the choice of appropriate chemotherapy.

*CTCs: Circulating Tumour Cells

ChemoSNiP

SNPs analysis of cytostatic drugs

Pharmacogenomics is the science that examines the inherited variations in genes that dictate drug response and explores the ways these variations can be used to predict whether a patient will have a good response to a drug, a bad response to a drug, or no response at all.

Single Nucleotide Polymorphisms (SNPs) are DNA sequence variations that can affect how human develop diseases or respond to chemicals, drugs and other agents.

ChemoSNiP includes the study of SNPs of genes that are included in metabolism/detoxification of cytostatic or targeted drugs which are used as cancer therapy.

Categories	Genes/Drugs	
Phase I Drug Metabolism	CYP2D6, CYP2C9, CYP1A2, CYP3A4	
	etc.	
Phase II Drug Metabolism	NAT, GSTP, TPMT etc.	
Pharmacodynamics	ABCB1, ABCG2 etc.	
	Oxaliplatin, Cisplatin,	
Alkylating Agents	Cyclophosphamide,	
	Melphalan, Busulfan, Dacarbazine etc.	
Topo I Inhibitors	Topotecan, Irinotecan etc.	
Topo II Inhibitors	Doxorubicin, Epirubicin, Etoposide etc.	
Antimetabolites	5-FU, Methotrexate, Gemcitabine,	
Spindle Poisons	Docetaxel, Paclitaxel, Vincristine etc.	

The polymorphisms are related with response to drugs, risk of toxicity, clearance, stability and activity of enzymes, reaction velocity, progression of free survival and overall survival.



SNPs analysis of cancer therapy and support drugs

PharmacoSNiP is based in the study of SNPs that related with the response or metabolism of broad spectrum of drug which are used in **cancer therapy** and **support**.

The test includes SNPs related with:

- Nonsteroidal anti-infammatory drugs (NSAIDs)
 - Antibiotics
 - Inflammation
 - Thrombosis risk
 - Diabetes risk
 - Osteoporosis risk
 - Heart Diseases risk, etc.

Polymorphisms which are associated with resistance, response, risk for gastrointestinal bleeding, clearance etc.

Indications of each test can be found at: http://www.rgcc-genlab.com/?tests

Where can I get these tests?

R.G.C.C. International GmbH



Research Genetic Cancer Centre

R.G.C.C. International GmbH is a leading company in analysis of circulating tumor cells as well as Cancer Stem Cells.

Through their analysis, is able to offer services in clinical field as well as in R&D in pharmaceutical industry.

By using the most advanced and innovative technologies of molecular and cellular biology, **R.G.C.C.**International GmbH manages to overpass several restrictions and difficulties that the analysis of CTCs (Circulating Tumor Cells) and CSCs (Circulating Stem Cells) involves.

Hence, through such an approach a massive amount of information and data has been generated in order to be used for identifying new drugable» targets as well as offering methods in clinical practice like new and precise assays, risk scale and classification of cancer patients.





R.G.C.C. International GmbH LOCATIONS OF GLOBAL NETWORK

RGCC INTERNATIONAL GmbH HEADQUARTERS

Baarerstrasse 95

Switzerland / Schweiz / Suisse

T.: +41 (0) 41 726 78 58,

F.: +41 (0) 41 726 78 59

E-mail: office@rgcc-international.com
Web address: www.rgcc-international.com

R.G.C.C. LTD. GREECE LAB FACILITIES

Industrial Area of Florina

GR 53100 - Florina, Greece

Tel.: +30 23850 41950, 41951, 41960, 41961

Fax.: +30 23850 41931

E-mail: office@rgcc-genlab.com Web address: www.rgcc-genlab.com

R.G.C.C. NORTH & WEST EUROPE Ltd BRANCH OFFICE IN NORTH & WEST EUROPE

Clifton, Bristol BS8 3LS, UK

T.: +44 (0) 117 317 1460, +44 (0) 779 505 6056

F.: +44 (0)117 973 3303

Email: info@rgcc-europe-nw.com

Web address: www.rgcc-europe-nw.com

R.G.C.C. USA LLC BRANCH OFFICE IN THE U.S.A., NORTH AMERICA CONTINENT & CANADA

3105 Main Street, Rowlett TX 75088, USA

T.: +1 214 299 9449

F.: +1 214 299 9923, +1 972 463 8243

Email: info@rgccusa.com

Web address: www.rgccusa.com

R.G.C.C. CENTRAL EUROPE GmbH BRANCH OFFICE IN CENTRAL EUROPE

Friedrichstrasse 88

10117 Berlin, Germany

T.: +49 30 81 00 37 51, +49 30 408173 281

F: +49 30 81 00 35 91

e-mail: office@rgcc-centraleurope.com

www.rgcc-centraleurope.com

MEDICAL GENETIC CANCER CENTER (M.G.C.C.) Ltd

BRANCH OFFICE IN CYPRUS, MIDDLE EAST & AFRICA

Achilleos 12A, Kato Lacatamia

2320 Nicosia, Cyprus

T.: +357 22 662869,

F.: +357 22 660104

E-mail: mgcccyprusltd@cytanet.com.cy

CANCER HOPE CLINIC REPRESENTATIVE IN ISRAEL

16 Hapaamonim st.

Raanana, 43391, Israel

T.: +972-3- 919 57-57

F.:+972-3-7256449

E-mail: info@cancerhope.co.il, office@cancerhope.co.il

Web address: www.cancerhope.co.il

DETOX CO Ltd REPRESENTATIVE IN JAPAN

Rokubancho State Building 103

5-18 Rokubancho Chiyodaku, Tokyo 102-0085, Japan

T..: 0081-(0)3-5876-4511

F.: 0081-(0)3-3222-5770 Email: info@detox.jp

Web address: www.detox.jp

NutriPATH REPRESENTATIVE IN OCEANIA

18A Harker St Burwood

Vic 3125 Australia

T.: +61 1300 688 522

F: +61 (03) 9880 2999

E-mail: info@nutripath.com.au

Web address: www.nutripath.com.au

1212 Limited

REPRESENTATIVE IN Hong Kong - China

25 Shum Wan Rd, Aberdeen, HONG KONG

T.: +852 25523108

F.: +852 30105801

Mob.: 00852 94887720

Email: drwaynemoran@1212.hk

www.1212.hk

Upmedical REPRESENTATIVE IN Brazil

Rua Cubatão, 395

<u>Pa</u>raiso,

04013-041 São Paulo/SP

Rrazil

T.: +55 11 4114 5962

Email: comercial@upmedical.com.br

Neeli Genetic Services REPRESENTATIVE IN INDIA

NSL Icon, 3rd Floor, Road Number 12, Plot No. 1-4,

Banjara Hills, Hyderabad, Telangana,

India 500034

T: (91) 789-321-1234

Email: office@neeligenetics.com

Website: www.neeligenetics.com

Where can i get these tests?